**35 PEEK CF Push-In
Suture Anchor with Suture**

Important Product Information
Produktinformationsblatt
Información importante sobre el Producto
Informations Importantes sur le Produit

Directions for Use
Gebrauchsanleitung
Instrucciones de uso
Mode d'emploi

30520 Rev B   0086
www.ParcusMedical.com
Mfd by: (941)755-7965 
Parcus Medical, LLC MDSS GmbH
6423 Parkland Dr. Schiffgraben 41
Sarasota, FL 34243 30175 Hanover, Germany

35 PEEK CF Push-In Suture Anchor with Suture

1. Indications

The Parcus 35 PEEK CF Push-In Suture Anchors with suture are indicated for attachment of soft tissue to bone. This product is intended for the following indications:

Shoulder Rotator Cuff Repair, Acromioclavicular Separation Repair, Bankart Lesion Repair, Biceps Tenodesis, Capsular Shift or Capsulolabral Reconstruction, Deltoid Repair, SLAP Lesion Repair.

Knee Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Posterior Oblique Ligament Repair, Extra Capsular Reconstruction, Iliotibial Band Tenodesis, Patellar Ligament and Tendon Avulsion Repair.

Foot/Ankle Lateral Stabilization, Medial Stabilization, Midfoot Reconstruction, Achilles Tendon Repair, Hallux Valgus Reconstruction, Metatarsal Ligament Repair.

Elbow Tennis Elbow Repair, Biceps Tendon Reattachment.

Hand/Wrist Scapholunate Ligament Reconstruction, Ulnar or Radial Collateral Ligament Reconstruction, TFCC.

2. Contraindications

- A. Any active infection.
- B. Blood supply limitations or other systemic conditions that may retard healing.
- C. Foreign body sensitivity, if suspected should be identified and precautions observed.
- D. Insufficient quality or quantity of bone. Suture anchor performance is directly related to the quality of bone into which the anchor is placed.
- E. Patient's inability or unwillingness to follow the surgeon's prescribed post-operative regimen.
- F. Any situation that would compromise the ability of the user to follow the directions for use or using the device for an indication other than those listed.

3. Adverse Effects

- A. Infection, both deep and superficial.
- B. Allergies and other reactions to device materials.
- C. Risks due to anesthesia.

4. Warnings

- A. This product should only be used by or on the order of a physician.
- B. The fixation provided by this device should be protected until healing is complete. Failure to follow the postoperative regimen prescribed by the surgeon could result in the failure of the device and compromised results.
- C. Anchor size selection should be made with care taking into consideration the quality of the bone in which the anchor is to be placed.
- D. Any decision to remove the device should take into consideration the potential risk of a second surgical procedure. Adequate postoperative management should be followed after implant removal.
- E. Polyethylene polyblend sutures, as used in this device, elicit a minimal inflammatory reaction in tissues, followed by gradual encapsulation of the suture by fibrous connective tissue. Polyethylene polyblend suture is not absorbed, nor is any significant change in tensile strength retention known to occur in vivo.
- F. The patient should be advised of the use and limitations of this device.
- G. Pre-operative planning and evaluation, surgical approaches and technique, and familiarity of the implant, including its instrumentation and limitations are necessary components in achieving a good surgical result.
- H. This device must never be reused. Reuse or re-sterilization may lead to changes in material characteristics such as deformation and material degradation which may compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection.
- I. This device must never be re-sterilized.
- J. Appropriate instrumentation should be used to implant this device.

5. Packaging and Labeling

- A. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature or has been altered in any way.

B. Please contact Parcus Medical Customer Service to report any package damage or alterations.

6. Material Specifications

The Parcus 35 PEEK CF Push-In Suture Anchors with suture is supplied attached to the high-strength, braided, polyethylene polyblend sutures either with or without stainless steel needles. The anchor material is polyetheretherkeytone (PEEK) Carbon Fiber Reinforced.

7. Sterilization

The Parcus 35 PEEK CF Push-In Suture Anchor with suture and driver assembly is supplied sterile. The contents are sterilized by EO gas. These products must never be re-sterilized.

8. Storage

Products must be stored in the original unopened package in a dry place and must not be used beyond the expiration date indicated on the package.

9. Directions for Use

A. Identify bone of sufficient quantity and quality into which the anchor is to be placed.

Note: In cases where bone quality is suspect and the 35 PEEK CF Push-In Suture Anchor with suture does not provide the desired fixation, the larger diameter threaded V-LoX™ PEEK CF Suture Anchors may provide an alternative solution.

- B. The Parcus 35 PEEK CF Push-In Suture Anchor with Suture is inserted into a socket created with one of two Parcus Drill Bits. The standard Drill Bit is 3.0mm. In less dense bone a 2.8mm Drill Bit should be used. Both the 2.8mm and 3.0mm Drill Bits are compatible with the Parcus Drill Guide. Arthroscopic cannulae, are typically used in conjunction with these devices. Always try to approach the targeted site as close to perpendicular as possible. If using the Drill Guide, use the obturator supplied with the Drill Guide to aid in its passage through tissue.
- C. Insert the Drill Bit into a drill chuck or pin driver to a depth that covers the laser etched line on the proximal end of the Drill Bit, as applicable. When using the Drill Guide, place the distal tip of the guide on the targeted bone and insert the Drill Bit through its lumen. Advance the Drill Bit into the bone approximately 16mm or until the circumferential laser mark on the distal tip of drill is flush with the surface of the bone. Remove the Drill Bit and Drill Guide, carefully noting the location of the prepared socket or leave the Drill Guide and Drill Bit in position until ready to insert the anchor.
- D. Aseptically open the package containing the 35 PEEK CF Push-In Suture Anchor with suture. Two strands of polyethylene polyblend high strength suture either with or without stainless steel needles are attached to the anchor. Note that each strand is a different color in order to aid in identification.
- E. If the Drill Bit was left in position as described in Step C, remove the Drill Bit before proceeding.
- F. Push the tip of the suture anchor into the prepared socket in line with the axis of the hole and using a suitable mallet, gently strike the driver handle, advancing the suture anchor into the prepared socket, until the circumferential laser mark on the distal end of the driver is flush with the bone.

Warning: Insertion of the anchor off axis from the hole may result in implant failure. Failure to insert the driver to the depth indicated by the laser mark may result in the anchor being left proud.

- G. If the 35 PEEK CF Push-In Suture Anchor with suture with needles is selected, release the suture from the back of the handle and disengage driver from anchor. The needles will now be loose and may be removed.
- H. Pass the sutures through the targeted tissue in the user preferred manner.
- I. Complete the tissue repair by tying knots of the surgeon preferred configuration and cutting the suture tails above the knots.

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