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Parcus PEEK and PEEK CF Suture Anchors

1. Indications:

The Parcus PEEK and PEEK CF Suture Anchors are indicated for attachment of soft tissue to bone. This product is intended for the following indications:

<u>Shoulder</u>	Rotator Cuff Repair, Acromioclavicular Separation Repair, Bankart Lesion Repair, Biceps Tenodesis, Capsular Shift or Capsulolabral Reconstruction, Deltoid Repair, SLAP Lesion Repair.
<u>Knee</u>	Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Posterior Oblique Ligament Repair, Extra Capsular Reconstruction, Iliotibial Band Tenodesis, Patellar Ligament and Tendon Avulsion Repair.
<u>Foot/Ankle</u>	Lateral Stabilization, Medial Stabilization, Midfoot Reconstruction, Achilles Tendon Repair, Hallux Valgus Reconstruction, Metatarsal Ligament Repair.
<u>Elbow</u>	Tennis Elbow Repair, Biceps Tendon Reattachment.
<u>Hand/Wrist</u>	Scapholunate Ligament Reconstruction, Ulnar or Radial Collateral Ligament Reconstruction, TFCC.

2. Contraindications:

- A. Any active infection.
- B. Blood supply limitations or other systemic conditions that may retard healing.
- C. Foreign body sensitivity, if suspected, should be identified and precautions observed.
- D. Insufficient quality or quantity of bone. Suture anchor performance is directly related to the quality of bone into which the anchor is placed.
- E. Patient's inability or unwillingness to follow the surgeon's prescribed post-operative regimen.
- F. Any situation that would compromise the ability of the user to follow the directions for use or using the device for an indication other than those listed.

3. Adverse Effects:

- A. Infection, both deep and superficial.
- B. Allergies and other reactions to device materials.
- C. Risks due to anesthesia.

4. Warnings:

- A. This product should only be used by or on the order of a physician.
- B. The fixation provided by this device should be protected until healing is complete. Failure to follow the postoperative regimen prescribed by the surgeon could result in the failure of the device and compromised results.
- C. Anchor size selection should be made with care taking into consideration the quality of the bone into which the anchor is to be placed. Osteopenic bone poses fixation challenges which may be addressed by larger diameter suture anchors.
- D. Any decision to remove the device should take into consideration the potential risk of a second surgical procedure. Adequate postoperative management should be followed after implant removal.
- E. Polyethylene polyblend sutures, as used in this device, elicit a minimal inflammatory reaction in tissues, followed by gradual encapsulation of the suture by fibrous connective tissue. Polyethylene polyblend suture is not absorbed, nor is any significant change in tensile strength retention known to occur in vivo.
- F. The patient should be advised of the use and limitations of this device.
- G. Pre-operative planning and evaluation, surgical approaches and technique, and familiarity of the implant, including its instrumentation and limitations are necessary components in achieving a good surgical result.
- H. This device must never be reused. Reuse or re-sterilization may lead to changes in material characteristics such as deformation and material degradation which may compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection.
- I. This device must never be re-sterilized.
- J. Appropriate instrumentation should be used to implant this device.

5. Packaging and Labeling:

- A. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature or has been altered in any way.
- B. Please contact Parcus Medical Customer Service to report any package damage or alterations.

6. Material Specifications:

PEEK and PEEK CF suture anchors are supplied with high-strength braided polyethylene polyblend sutures with or without stainless steel needles. The suture anchors are manufactured from Polyetheretherkeytone with or without carbon fiber.

7. Sterilization:

The PEEK and PEEK CF Suture Anchor and driver assemblies are supplied sterile. These products must never be re-sterilized.

8. Storage:

Products must be stored in the original unopened package in a dry place and must not be used beyond the expiration date indicated on the package.

9. Directions for Use

- A. Identify bone of sufficient quantity and quality into which the anchor is to be placed.

Note: When bone quality is compromised, smaller diameter suture anchors may not provide adequate fixation. The use of larger diameter suture anchors may provide an alternative solution.

- B. Select the desired configuration (either with or without needles) and aseptically open the package containing the suture anchor.
- C. Taking care to create an entry to the desired placement site as close to perpendicular as possible, insert the appropriate hole preparation instrument per associated instrumentation DFU. The instrumentation associated with each size anchor can be found in the table below.

Note: Smaller diameter hole preparation instruments are used to increase fixation in less dense bone.

Anchor	Instrumentation
4.5mm Twist	3.2mm Awl or Drill
	4.0mm Drill
5.5mm Twist or V-LoX PEEK CF	Parcus 5.5 – 6.5mm Punch/Tap
6.5mm Twist or V-LoX PEEK CF	Parcus 5.5 – 6.5mm Punch/Tap

- D. Place the PEEK or PEEK CF Suture Anchor into the prepared site, as described above, and turn the driver handle clockwise until the circumferential laser etched line on the distal end of the driver shaft is flush with the surrounding bone. If present the longitudinal lines on the driver shaft are used to indicate the orientation of the suture passing around the post within the anchor. Placing the longitudinal lines in a plane perpendicular to the targeted tissue offers the best orientation to allow the suture to slide through the anchor.
- E. If a suture anchor with needles is selected, release the suture loops from the back of the handle and disengage driver from anchor. The needles will now be loose and may be removed.
- F. Two or three strands of polyethylene polyblend high strength suture either with or without needles are attached to the anchor. Note that each strand is a different color in order to aid in identification.
- G. Pass the sutures through the targeted tissue in the user preferred manner.
- H. Complete the tissue repair by tying knots of the surgeon preferred configuration and cutting the suture tails above the knots.

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