**Parcus Drill Guide  
for Suture Anchor Placement**

**Important Product Information**

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## Parcus Drill Guide for Suture Anchor Placement

### 1. Indications:

Parcus Drill Guides are provided in various configurations to assist in the placement of suture anchors and/or drill bits.

### 2. Warnings:

- A. This product is intended for use by or on the order of a physician.
- B. Parcus Drill Guides are supplied unsterile as reusable devices.
- C. Never use a bent or damaged drill guide or obturator.
- D. This product is not intended for use as an implant.
- E. Loading this instrument or any instrument used with this device in an off axis manner may result in breakage or bending of the instrument, drill guide or obturator.

### 3. Material:

This device is comprised of a Stainless Steel Alloy shaft and an Aluminum handle. The materials used in the manufacture of this device that are intended to be placed inside the body are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

### 4. Inspection:

- A. Inspect the device for damage at all stages of handling.
- B. If damage is detected, consult the manufacturer for guidance.

### 5. Cleaning:

- A. Immediate rinsing and cleaning after use with an enzymatic detergent will effectively remove and prevent drying of adherent blood, tissue, etc.
- B. Scrub device with a soft brush, paying special attention to areas where debris might accumulate. Always avoid any harsh materials that can scratch or mar the surface of the device.
- C. Rinse the device thoroughly with water following the cleaning process.
- D. The device should be placed in an ultrasonic cleaning unit and processed according to the unit's directions.
- E. The device should be rinsed thoroughly with water following the ultrasonic process.

### 6. Sterilization:

This device is not provided sterile. The following table provides recommended minimum sterilization parameters that have been validated by Parcus Medical to provide a  $10^{-6}$  sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Gravity	132°C / 270°F	15 minutes	30 minutes <sup>1</sup>
Pre-vacuum	132°C / 270°F	4 minutes	

<sup>1</sup> Drying times vary according to load size and should be increased for larger loads.

### 7. Directions for use:

- A. Place the Drill Guide with obturator, as applicable, on the targeted bone through an open incision, percutaneously, or through an arthroscopic cannula.
- B. Remove obturator, as applicable, and advance the required instrumentation through the guide to the desired depth.
- C. Remove the instrumentation from the guide.
- D. Remove the Drill Guide, as necessary.

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