**Parcus Braid™
Non-Absorbable Surgical Suture**

Important Product Information

Produktinformationsblatt
Información Importante sobre el Producto
Informations Importantes sur le Produit
Önemli Ürün Bilgileri
Informazioni Importanti sul Prodotto
Informações Importantes sobre o Produto

Directions for Use

Gebrauchsanleitung
Instrucciones de uso
Mode d'emploi
Kullanım talimatları
Istruzioni per l'uso
Instruções para o Uso

30838 Rev 4

Mfd by:
Parcus Medical, LLC
6423 Parkland Dr.
Sarasota, FL 34243

www.ParcusMedical.com
(941)755-7965

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Parcus Braid™
Non-Absorbable Surgical Suture

DESCRIPTION: Parcus Braid sutures are non-absorbable, sterile, surgical sutures composed of ultra high molecular weight polyethylene (UHMWPE). The monofilament fibers are uncoated and the braided sutures are available uncoated or coated uniformly with polybutylene adipate or silicone to increase surface lubricity, thereby enhancing the handling characteristics, ease of passage through tissue, and knot run-down properties for security. Parcus Braid sutures are inert and elicit only minimal local tissue reaction. The suture braid is available undyed (white) or with trace filaments of black nylon or blue polypropylene suture for color. Parcus Braid sutures meet all requirements established by the United States Pharmacopeia (USP) and the European Pharmacopeia (EP) for non-absorbable surgical sutures except where indicated on the individual unit package.

INDICATIONS: Parcus Braid sutures are indicated for use in general soft tissue approximation and/or ligation and the use of allograft tissues for orthopedic surgeries. Parcus Braid sutures are intended for one-time use only, and are not to be re-sterilized.

ACTIONS: Parcus Braid sutures elicit a minimal acute inflammatory reaction in tissues, followed by gradual encapsulation of the suture by fibrous connective tissue. Parcus Braid sutures are not absorbed, nor is there any significant change in tensile strength retention known to occur in vivo.

CONTRAINDICATIONS: None known.

WARNINGS: As with any foreign body, prolonged contact of any suture with salt solutions may result in calculus formation. Do not resterilize. Discard opened or unused sutures. Users should be familiar with surgical procedures and techniques involving sutures before employing Parcus Braid non-absorbable sutures for wound closure, as the risk of wound dehiscence may vary with the site of application and the suture material used. Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

PRECAUTIONS: In handling this or any other suture material, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders. Adequate knot security requires the accepted surgical technique of flat, square ties, with additional throws as warranted by surgical circumstance and the experience of the surgeon. The use of additional throws may be particularly appropriate when knotting monofilaments.

ADVERSE REACTIONS: Adverse effects associated with the use of this device include: wound dehiscence, calculi formation when prolonged contact with salt solutions occur, enhanced bacterial infectivity, minimal acute inflammatory reaction, and transitory local irritation.

HOW SUPPLIED: Parcus Braid sutures are available undyed (white) or with trace filament of black nylon or blue polypropylene braided into the suture for color. Parcus Braid is available coated or uncoated, in sizes 7 (9 metric) through 6-0 (0.7 Metric). Parcus Suture Tape is a flat braid suture offered in widths from 1 mm to 5 mm. The sutures are supplied sterile, in pre-cut lengths and ligating reels, non-needled or affixed to needles using both permanent and removable needle attachment techniques. Parcus Braid sutures are available in 0.5, one, two, and three dozen box quantities.

CAUTION: Federal (USA) law restricts this device to sale and use by, or on the order of, a physician.

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