Parcus SLiK Anchor

Important Product Information

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Parcus SLiK Anchor (English)

1. Indications

The Parcus SLiK Anchors are indicated for attachment of soft tissue to bone. This product is intended for the following indications:

Shoulder Rotator Cuff Repair, Acromioclavicular Separation Repair, Bankart Lesion Repair, Biceps Tenodesis, Capsular Shift or Capsulolabral Reconstruction, Deltoid Repair, SLAP Lesion Repair.

Knee Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Posterior Oblique Ligament Repair, Extra Capsular Reconstruction, Iliotibial Band Tenodesis, Patellar Ligament and Tendon Avulsion Repair.

Foot/Ankle Lateral Stabilization, Medial Stabilization, Midfoot Reconstruction, Achilles Tendon Repair, Hallux Valgus Reconstruction, Metatarsal Ligament Repair.

Elbow Tennis Elbow Repair, Biceps Tendon Reattachment.

Hand/Wrist Scapholunate Ligament Reconstruction, Ulnar or Radial Collateral Ligament Reconstruction, TFCC.

2. Contraindications

- A. Any active infection.
- B. Blood supply limitations or other systemic conditions that may retard healing.
- C. Foreign body sensitivity, if suspected should be identified and precautions observed.
- D. Insufficient quality or quantity of bone. Suture anchor performance is directly related to the quality of bone into which the anchor is placed.
- E. Patient's inability or unwillingness to follow the surgeon's prescribed post-operative regimen.
- F. Any situation that would compromise the ability of the user to follow the directions for use or using the device for an indication other than those listed.

3. Adverse Effects

- A. Infection, both deep and superficial.
- B. Allergies and other reactions to device materials.
- C. Risks due to anesthesia.

4. Warnings

- A. This product should only be used by or on the order of a physician.
- B. The fixation provided by this device should be protected until healing is complete. Failure to follow the postoperative regimen prescribed by the surgeon could result in the failure of the device and compromised results.
- C. Anchor size selection should be made with care taking into consideration the quality of the bone in which the anchor is to be placed.
- D. Any decision to remove the device should take into consideration the potential risk of a second surgical procedure. Adequate postoperative management should be followed after implant removal.
- E. The patient should be advised of the use and limitations of this device.

- F. Pre-operative planning and evaluation, surgical approaches and technique, and familiarity of the implant, including its instrumentation and limitations are necessary components in achieving a good surgical result.
- G. This device must never be reused. Reuse or re-sterilization may lead to changes in material characteristics such as deformation and material degradation which may compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection.
- H. This device must never be re-sterilized.
- I. Appropriate instrumentation should be used to implant this device.

5. Packaging and Labeling

- A. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature or has been altered in any way.
- B. Please contact Parcus Medical Customer Service to report any package damage or alterations.

6. Material Specifications

The Parcus SLiK Anchors are manufactured from PEEK (Polyetheretherketone) and/or PEEK CF (Polyetheretherketone carbon fiber reinforced). These devices are MR Safe.

7. Sterilization

The Parcus SLiK Anchors are supplied in sterile packaging. The contents are sterilized by EO gas. These products must never be re-sterilized.

8. Storage

Products must be stored in the original unopened package in a dry place and must not be used beyond the expiration date indicated on the package.

9. Directions for Use

- A. Identify bone of sufficient quantity and quality into which the anchor is to be placed.
- B. The Parcus SLiK Anchor is inserted into a socket created with a Parcus Drill Bit or Parcus Awl. See the table below for drill bit and awl selection.

Instrumentation	Used With
4.0mm Parcus Drill Bit (Reference # 11000U)	4.75x15 mm SLiK Anchor 5.5x15 mm SLiK Anchor 6.25x15 mm SLiK Anchor
4.1mm Parcus Awl (Reference # 11092)	4.75x15 mm SLiK Anchor 5.5x15 mm SLiK Anchor 6.25x15 mm SLiK Anchor
4.8mm Parcus Awl (Reference # 11093)	5.5x15 mm SLiK Anchor 6.25x15 mm SLiK Anchor

- C. Arthroscopic cannulae are typically used in conjunction with these devices. Always try to approach the targeted site as close to perpendicular as possible. If using a Drill Guide, use the obturator supplied with the Drill Guide to aid in its passage through tissue.
- D. Place sutures (not supplied with anchor) through the targeted tissue using the user's preferred method.
- E. Use the Awl or Drill Bit to create a socket into which the anchor is to be placed. When using an Awl, identify the site for the placement of the anchor and use a mallet to drive the awl to the depth

indicated by the the circumferential line near its distal tip. This mark must be inserted flush with the surface of the targeted bone. When using a Drill Bit, place the bit into a drill chuck or pin driver to a depth that covers the laser etched line or collar stop on the proximal end of the Drill Bit, as applicable. Insert the drill until the circumferential line, located near the distal tip of the drill, is flush with the surface of the targeted bone. When using the Drill Guide, place the distal tip of the guide on the targeted bone and insert the Drill Bit through its lumen. Advance the Drill Bit into the bone approximately 25mm or until the circumferential laser mark on the distal tip of the Drill Bit is flush with the bone. If the drill bit features a collar stop, advance the drill bit into the bone until the collar stop contacts the proximal end of the Guide. Remove the Drill Bit and Drill Guide, carefully noting the location of the prepared socket or leave the Drill Guide and Drill Bit/Awl in position until ready to insert the anchor.

- F. Aseptically open the package containing the Parcus SLiK Anchor. Insure that the anchor body is fully engaged on the shoulder of the driver by applying slight force on the anchor body in a proximal direction. The tails of the sutures that were previously placed through the target tissue are then passed through the Suture Guide located at the distal tip of the SLiK Anchor Driver. The Parcus SLiK Anchor will accommodate up to two (2) strands of 2mm wide Suture Tape and one (1) strand of #2 USP suture. In arthroscopic applications, draw the sutures up through the same cannula through which the Drill Guide will be used.
- G. Place the distal tip of the Suture Guide into the the previously created bone socket until the anchor body makes contact with the surface of the bone. Apply force to the suture tails until the target tissue has achieved the desired tension.
- H. Depress the button on the driver handle and screw the anchor body into the bone by applying force and rotating in a clock-wise direction. Insert until the circumferential laser mark on the distal tip of the driver is flush with the surface of the bone.
Warning: Insertion of the anchor off axis from the hole may result in implant failure. Failure to insert the driver to the depth indicated by the laser mark may result in the anchor being left proud.
- I. Remove the driver, internal driver shaft, and Suture Guide by pulling the driver straight out. It is noted that the Suture Guide has a necked down section that will allow the passed suture to break through and permit the Suture Guide to be pulled through the ID of the anchor body.
- J. Cut the ends of the remaining tails of the suture.

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