## Parcus Graft Sizer

### Important Product Information

Produktinformationsblatt

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MDSS GmbH

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30175 Hanover, Germany

**Parcus Graft Sizer**  
**Product Information Sheet and Directions for Use:**

**1. Indications:**

The Parcus Graft Sizers are used to assist in determining the diameter of allograft and autograft tissue used for reconstructing the Anterior and Posterior Cruciate Ligaments.

**2. Warnings:**

- A. This product is intended for use by or on the order of a physician.
- B. No component of this instrument is intended to be left in the body.

**3. Material:**

This product is manufactured from aluminum or a stainless steel alloy.

**4. Inspection:**

- A. Inspect the device for damage at all stages of handling.
- B. If damage is detected, consult the manufacturer for guidance.

**5. Cleaning:**

- A. Immediate rinsing and cleaning after use with an enzymatic detergent will effectively remove and prevent drying of adherent blood, tissue, etc.
- B. Scrub device with a soft brush, paying special attention to the cannulation and other areas where debris might accumulate. Always avoid any harsh materials that can scratch or mar the surface of the device.
- C. Rinse the device thoroughly with water following the cleaning process.

**6. Sterilization:**

This product is a non-sterile item that must be cleaned and sterilized prior to use. The following table provides recommended minimum sterilization parameters that have been validated by Parcus Medical to provide a  $10^{-6}$  sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Gravity	132°C / 270°F	15 minutes	30 minutes <sup>1</sup>
Pre-vacuum	132°C / 270°F	4 minutes	

<sup>1</sup> Drying times vary according to load size and should be increased for larger loads

**7. Directions for use:**

- A. Using the sutures attached to a properly prepared graft, draw the graft into the sizing sleeve. The sleeve is used as Go-No Go gauge.
- B. The diameter of the graft will be determined when it passes snugly through the sleeve.

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