## Parcus SLiK Fix Screw

### Important Product Information

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## Parcus SLiK Fix Screw (English)

### 1. Indications

The Parcus SLiK Fix Screws are intended for use in soft tissue to bone reattachment procedures. Specific indications are the following:

Shoulder Rotator Cuff Repair, Acromioclavicular Separation Repair, Bankart Lesion Repair, Biceps Tenodesis, Capsular Shift or Capsulolabral Reconstruction, Deltoid Repair, SLAP Lesion Repair.

Knee Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Posterior Oblique Ligament Repair, Extra Capsular Reconstruction, Iliotibial Band Tenodesis, Patellar Ligament and Tendon Avulsion Repair.

Foot/Ankle Lateral Stabilization, Medial Stabilization, Midfoot Reconstruction, Achilles Tendon Repair, Hallux Valgus Reconstruction, Metatarsal Ligament Repair.

Elbow Tennis Elbow Repair, Biceps Tendon Reattachment.

Hand/Wrist Scapholunate Ligament Reconstruction, Ulnar or Radial Collateral Ligament Reconstruction, TFCC.

### 2. Contraindications

- A. Any active infection.
- B. Blood supply limitations or other systemic conditions that may retard healing.
- C. Foreign body sensitivity, if suspected should be identified and precautions observed.
- D. Insufficient quality or quantity of bone. Suture anchor performance is directly related to the quality of bone into which the anchor is placed.
- E. Patient's inability or unwillingness to follow the surgeon's prescribed post-operative regimen.
- F. Any situation that would compromise the ability of the user to follow the instructions for use or using the device for an indication other than those listed.

### 3. Adverse Effects

- A. Infection, both deep and superficial.
- B. Allergies and other reactions to device materials.
- C. Risks due to anesthesia.

### 4. Warnings

- A. This product should only be used by or on the order of a physician.
- B. The fixation provided by this device should be protected until healing is complete. Failure to follow the postoperative regimen prescribed by the surgeon could result in the failure of the device and compromised results.
- C. Screw size selection should be made with care taking into consideration the quality of the bone in which the anchor is to be placed.
- D. Any decision to remove the device should take into consideration the potential risk of a second surgical procedure. Adequate postoperative management should be followed after implant removal.
- E. The patient should be advised of the use and limitations of this device.

- F. Pre-operative planning and evaluation, surgical approaches and technique, and familiarity of the implant, including its instrumentation and limitations are necessary components in achieving a good surgical result.
- G. This device must never be reused. Reuse or re-sterilization may lead to changes in material characteristics such as deformation and material degradation which may compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection.
- H. This device must never be re-sterilized.
- I. Appropriate instrumentation should be used to implant this device.

## **5. Packaging and Labeling**

- A. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature or has been altered in any way.
- B. Please contact Parcus Medical Customer Service to report any package damage or alterations.

## **6. Material Specifications**

The Parcus SLiK Fix Screws are manufactured from PEEK (Polyetheretherketone) and/or PEEK CF (Polyetheretherketone carbon fiber reinforced). The tissue protection washer on the insertion and tensioning rod is manufactured from PEEK (Polyetheretherketone). These devices are MR Safe.

## **7. Sterilization**

The Parcus SLiK Fix Screws are supplied in sterile packaging. The contents are sterilized by EO gas. These products must never be re-sterilized.

## **8. Storage**

Products must be stored in the original unopened package in a dry place and must not be used beyond the expiration date indicated on the package.

## **9. Instructions for Use**

- A. Identify bone of sufficient quantity and quality into which the screw is to be placed.
- B. The Parcus SLiK Fix Screw is inserted into a bone socket using the user preferred technique.
- C. Arthroscopic cannulae are typically used in conjunction with these devices. Always try to approach the targeted site as close to perpendicular as possible. If using a Drill Guide, use the obturator supplied with the Drill Guide to aid in its passage through tissue.
- D. When using the SLiK Fix for suture fixation, place sutures (not supplied with anchor) through the targeted tissue using the user's preferred method.
- E. Select a screw of the appropriate diameter to achieve the desired fixation. For suture fixation, the screw should be 1mm larger than the diameter of the hole. For graft fixation, the screw should be 1mm smaller than the diameter of the hole.
- F. Aseptically open the package containing the Parcus SLiK Fix Screw. Install the anchor body into full engagement with the driver by applying slight force on the anchor body in a proximal direction.
- G. Use the trocar tip on the distal end of the Insertion and Tensioning Rod to capture the graft or suture in the desired location.
- H. Place the distal tip of the Insertion and Tensioning Rod into the bone socket until the screw makes contact with the surface of the bone. Apply tension to the suture tails as necessary.
- I. Depress the button on the handle and insert the anchor body into the bone by applying force in the axial direction and rotating clock-wise. Insert until the circumferential laser mark on the distal end of the driver is flush with the surface of the bone.

**Warning: Insertion of the anchor off axis from the hole may result in implant failure. Failure to insert the driver to the depth indicated by the laser mark may result in the anchor being left proud.**

- J. Remove the Driver with the Insertion and Tensioning Rod by pulling the driver straight out. The PEEK Tissue Protection Washer on the distal tip will be pressed off of the internal rod and remain implanted in the bone socket.
- K. When necessary, cut off the tails of the suture using the preferred suture cutting device.

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