**Directions for Use
Parcus Chondral Picks**

- 1. Indications:**
The Parcus Chondral Picks are designed to create a socket/defect in bone.
- 2. Warnings:**
 - A. This product is intended for use by or on the order of a physician.
 - B. Do not use this product if it is bent or there is visible damage to the shaft or drive section.
 - C. This device is not intended to be used as an implant.
- 3. Material:**
This device is manufactured from a Stainless Steel Alloy and may include an aluminum alloy handle. The materials used in the manufacture of this device are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.
- 4. Inspection:**
 - A. Inspect the device for damage at all stages of handling.
 - B. If damage is detected, consult the manufacturer for guidance.
- 5. Cleaning:**
 - A. Immediate rinsing and cleaning after use with an enzymatic detergent will effectively remove and prevent drying of adherent blood, tissue, etc.
 - B. Scrub device with a soft brush, paying special attention to areas where debris might accumulate. Always avoid any harsh materials that can scratch or mar the surface of the device.
 - C. Rinse the device thoroughly with water following the cleaning process.
- 6. Sterilization:**
This device is not provided sterile. The following table provides recommended minimum sterilization parameters that have been validated by Parcus Medical to provide a 10⁻⁶ sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Gravity	132°C / 270°F	15 minutes	30 minutes ¹
Pre-vacuum	132°C / 270°F	4 minutes	

¹ Drying times vary according to load size and should be increased for larger loads.

- 7. Directions for use:**
 - A. Place the tip of the Chondral Pick on the surface of the bone at the point at which the socket/defect is to be created.
 - B. As necessary use a suitable mallet to insert the Chondral Pick into the bone to the desired depth. The circumferential laser line indicates a distance of 3.5 mm from the tip.
 - C. Remove the Chondral Pick

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