Twist Knotless & V-LoX Knotless Anchors

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Twist Knotless & V-LoX Knotless Anchors (English)

1. Indications:

The Parcus Twist Knotless & V-LoX Knotless Anchors are indicated for attachment of soft tissue to bone. These products are intended for the following indications:

- Shoulder: Rotator Cuff Repair, Acromioclavicular Separation Repair, Bankart Lesion Repair, Biceps Tenodesis, Capsular Shift or Capsulolabral Reconstruction, Deltoid Repair, SLAP Lesion Repair.
- Knee: Medial Collateral Ligament Repair, Lateral Collateral Ligament Repair, Posterior Oblique Ligament Repair, Extra Capsular Reconstruction, Iliotibial Band Tenodesis, Patellar Ligament and Tendon Avulsion Repair.
- Foot / Ankle: Lateral Stabilization, Medial Stabilization, Midfoot Reconstruction, Achilles Tendon Repair, Hallux Valgus Reconstruction, Metatarsal Ligament Repair.
- Elbow: Tennis Elbow Repair, Biceps Tendon Reattachment.
- Hand / Wrist: Scapholunate Ligament Reconstruction, Ulnar or Radial Collateral Ligament Reconstruction, TFCC.

2. Contraindications:

- A. Any active infection.
- B. Blood supply limitations or other systemic conditions that may retard healing.
- C. Foreign body sensitivity, if suspected should be identified and precautions observed.
- D. Insufficient quality or quantity of bone. Suture anchor performance is directly related to the quality of bone into which the anchor is placed.
- E. Patient's inability or unwillingness to follow the surgeon's prescribed post-operative regimen.
- F. Any situation that would compromise the ability of the user to follow the instructions for use or using the device for an indication other than those listed.

3. Adverse Effects:

- A. Infection, both deep and superficial.
- B. Allergies and other reactions to device materials.
- C. Risks due to anesthesia.

4. Warnings:

- A. This product should only be used by or on the order of a physician.
- B. The fixation provided by this device should be protected until healing is complete. Failure to follow the postoperative regimen prescribed by the surgeon could result in the failure of the device and compromised results.
- C. Anchor size selection should be made with care taking into consideration the quality of the bone in which the anchor is to be placed.
- D. Any decision to remove the device should take into consideration the potential risk of a second surgical procedure. Adequate postoperative management should be followed after implant removal.
- E. The patient should be advised of the use and limitations of this device.
- F. Pre-operative planning and evaluation, surgical approaches and technique, and familiarity of the implant, including its instrumentation and limitations are necessary components in achieving a good surgical result.
- G. This device must never be reused. Reuse or re-sterilization may lead to changes in material characteristics such as deformation and material degradation which may compromise device

performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection.

H. This device must never be re-sterilized.

I. Appropriate instrumentation should be used to implant this device.

5. Packaging and Labeling:

A. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature or has been altered in any way.

B. Please contact Parcus Medical Customer Service to report any package damage or alterations.

6. Material Specifications:

The Parcus Twist Knotless is manufactured from PEEK (Polyetheretherketone). The Parcus V-LoX Knotless Anchor is manufactured from PEEK CF (Polyetheretherketone Carbon Fiber reinforced). These devices are MR Safe.

7. Sterilization:

The Parcus Twist Knotless & V-LoX Knotless Anchors are supplied in sterile packaging. The contents are sterilized by EO gas. These products must never be re-sterilized.

8. Storage:

Products must be stored in the original unopened package in a dry place and must not be used beyond the expiration date indicated on the package.

9. Instructions for Use:

A. Identify bone of sufficient quantity and quality into which the anchor is to be placed.

B. The Parcus Twist Knotless & V-LoX Knotless Anchors are inserted into a socket created with a Parcus Drill Bit or Parcus Awl. See the table below for drill bit and awl selection.

Instrumentation Used with Twist Knotless & V-LoX Knotless Suture Anchors

Note: Drill and Awl diameter are dependent upon bone density. Generally, less dense or osteopenic bone will use the smallest drill or awl diameter listed.

Instrument	Reference #	Anchor Diameter
4.0 mm Parcus Drill Bit	11000U	4.75, 5.5, 6.25 mm
Drill Guide 100 mm with Beaver Tail Handle	11241	N/A
Tapered Awl	11172	4.75, 5.5, 6.25 mm

C. Arthroscopic cannulae are typically used in conjunction with these devices. Always try to approach the targeted site as close to perpendicular as possible.

D. Place sutures (not supplied with anchor) through the targeted tissue using the user's preferred method.

E. Use the Awl or Drill Bit to create a socket into which the anchor is to be placed. When using an Awl, identify the site for the placement of the anchor and use a mallet to drive the awl to the depth indicated by one of the circumferential lines near its distal tip. **(Note: The distal line is used in less dense bone. The proximal line is used in more dense bone.)** This mark must be inserted flush with the surface of the targeted bone. When using the drill bit without the drill guide, advance the drill so that the most distal circumferential line is flush with the surface of the targeted bone. When using the drill bit with the drill guide, place the distal tip of the guide on the targeted bone and insert the Drill Bit through its lumen. Advance the drill bit until the most proximal laser line, (adjacent to the stop collar) is flush with the proximal end of the guide. Remove the drill bit and drill guide, carefully noting the location of the prepared socket. The drill guide and drill bit/awl may be left in position to maintain visualization of the socket until the anchor is ready for insertion.

- F. Aseptically open the package containing the Parcus Twist Knotless or V-LoX Knotless Anchor. Ensure that the anchor body is fully engaged on the shoulder of the driver by applying slight force on the anchor body in a proximal direction. The tails of the sutures that were previously placed through the target tissue are then passed through the eye of the wire suture threader located at the distal tip of the Twist Knotless or V-LoX Knotless Anchor Driver (the Parcus Anchor will accommodate up to four (4) strands of Suture Tape or four (4) strands of #2 USP suture).
- G. Pull on the handle of the suture threader protruding from the proximal end of the driver to pull the suture or tape up through the cannulated shaft. Pull the desired amount of tension on the suture and then insert the tip of the anchor body into the previously created bone socket. Turn the driver in a clockwise direction, while releasing tension on the sutures, to advance the anchor body until the circumferential laser mark on the distal tip of the driver is flush with the surface of the surrounding bone.

Warning: Insertion of the anchor off axis from that of the socket may result in implant failure. Failure to insert the driver to the depth indicated by the laser mark may result in the anchor being left proud.

- H. Remove the driver by pulling the driver straight out.
- I. Cut the ends of the remaining tails of the suture.

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