Parcus GFS Naked

Important Product Information
Produktinformationsblatt
Información Importante sobre el Producto
Informations Importantes sur le Produit

Instructions for Use
Gebrauchsanleitung
Instrucciones de uso
Mode d'emploi

31152 Rev 2
Eff. 03/20

 **PARCUS**
www.ParcusMedical.com
(941)755-7965

 Mfd by:
Parcus Medical, LLC
6423 Parkland Dr.
Sarasota, FL 34243

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Parcus Graft Fixation System Naked Product Family (English)

1. Indications:

The Parcus GFS Naked devices are indicated for use in the fixation of ligaments and tendons in patients requiring ligament or tendon repair.

2. Contraindications:

- A. Any active infection.
- B. Blood supply limitations or other systemic conditions that may retard healing.
- C. Foreign body sensitivity, if suspected, should be identified and precautions observed.
- D. Insufficient quality or quantity of bone.
- E. Patient's inability or unwillingness to follow surgeon's prescribed post-operative regimen.
- F. Any situation that would compromise the ability of the user to follow the directions for use or using the device for an indication other than those listed.
- G. The GFS Naked devices are intended to be used in conjunction with a suture product such as #2 or #5 Parcus Braid Suture which is ultra high-strength, braided, polyethylene polyblend sutures. Use of this device with other types or sizes of suture products may lower the ultimate strength of the device which could compromise the overall fixation. Implant size and suture selection is at the discretion of the surgeon and should be based on the expected forces that the fixation is intended to withstand.

3. Adverse Effects:

- A. Infection, both deep and superficial.
- B. Allergies and other reactions to device materials.
- C. Risks due to anesthesia.

4. Warnings:

- A. Caution: Federal law restricts this device to sale by or on the order of a physician.
- B. The fixation provided by this device should be protected until healing is complete. Failure to follow the postoperative regimen prescribed by the surgeon could result in the failure of the device and compromised results.
- C. Size selection of the implant and suture should be made with care taking into consideration the quality of the bone on which the implant will rest and the anticipated forces that the fixation must resist.
- D. Any decision to remove the device should take into consideration the potential risk of a second surgical procedure. An adequate postoperative management plan should be implemented after implant removal.
- E. Pre-operative planning and evaluation, surgical approaches and technique, and familiarity of the implant, including its instrumentation and limitations are necessary components in achieving a good surgical result.
- F. This device must never be reused. Reuse or re-sterilization may lead to changes in material characteristics such as deformation and material degradation which may compromise device performance. Reprocessing of single use devices can also cause cross-contamination leading to patient infection.
- G. This device must never be re-sterilized.
- H. Appropriate instrumentation should be used to implant this device.

- I. This device has not been evaluated for safety and compatibility in the MR environment. This device has not been tested for heating or migration in the MR environment. Use of MR technology in the presence of devices of this nature may cause magnetically induced displacement forces and torques, radio frequency heating and image artifacts. Standard MRI screening guidelines for post-operative patients should be followed.

5. Packaging and Labeling:

- A. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature or has been altered in any way.
- B. Please contact Parcus Medical Customer Service to report any package damage or alterations.

6. Material Specifications:

The Parcus GFS Naked Product Family is comprised of Ti-6Al-4V ELI (ASTM F136).

7. Sterilization:

This Parcus GFS Naked Product Family is supplied sterile.

8. Storage:

Products must be stored in the original unopened package in a dry place and must not be used beyond the expiration date indicated on the package.

9. Directions for Use (Soft Tissue Graft)

- A. A bone tunnel, appropriately sized to accommodate the graft and applicable GFS Naked, is necessary for use of the GFS Naked device. Based on the implant chosen and size of graft used, a stepped tunnel may be required. Please refer to product packaging for recommended diameter for distal cortex.
- B. Using an appropriately sized Drill Tip Guide Pin or Drill, create bone tunnel through to the outer cortex of the bone. If tunnel is within the specified size of the implant per package labeling and large enough to accommodate the graft, no further drilling or reaming is required. However, if tunnel is not large enough, use of a cannulated drill or reamer is recommended. If utilizing a stepped tunnel, use the Parcus GFS Depth Gauge to measure the tunnel in order to ensure that the lateral cortex of the bone is not breached as that may compromise the strength of the fixation.

Note: If the lateral cortex of the bone is breached, use of a larger GFS Naked device may be required.

- C. Once the graft has been accurately measured using a Parcus Graft Sizer or Graft Sizing Block and the tunnel has been prepared, the GFS Naked may be assembled to the graft.
- D. Whipstitch the graft in the surgeon preferred manner, leaving suture tails long enough to be placed through the two center holes of the GFS Naked fixation device. Aseptically open the appropriately sized Parcus GFS Naked device and pass suture tails through center holes of implant.

Note: The GFS Nano Naked can only accommodate suture up to #2 in size. However, all larger sizes of GFS Naked can accommodate suture up to #5 in size.

- E. Push the GFS Naked through the tunnel using an appropriate instrument and then draw the graft into position by pulling on the sutures lying along the length of the graft or pull the device through the tunnel by adding accessory sutures, one each to either end of the device. Pass the sutures

through the tunnel and then pull the device to a point where the device exits from the cortex distal to the graft, taking care not to pull the device into the surrounding soft tissue.

- F. Complete the fixation by passing the sutures back through the graft and tying knots or by using adjunctive fixation, e.g. placing the sutures around a bicortical screw and washer.

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