**Parcus Drill Tip Guide Pin  
With or Without Suture Eyelet**

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## Parcus Drill Tip Guide Pin With or Without Suture Eyelet (English)

### 1. Indications:

This pin is used to aid in the guidance of cannulated drills and reamers. The Drill Tip Guide Pin with Suture Eyelet can be effectively used to assist in the passing of suture.

### 2. Warnings:

- A. Caution: Federal Law restricts this device to sale by or on the order of a physician.
- B. This product is designed for single use only.
- C. Never use a bent pin or a pin that has surface defects. Pins that are reused can fatigue and break.
- D. This product is not intended for use as an implant.
- E. Use of drill or reamer off-axis from guide pin can damage pin and result in product failure.

### 3. Material:

This device is manufactured from a Stainless-Steel Alloy. The materials used in the manufacture of this device that are intended to be placed inside the body are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

### 4. Sterilization:

This device is not provided sterile. The following table provides recommended minimum sterilization parameters that have been validated by Parcus Medical to provide a  $10^{-6}$  sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Gravity	132°C / 270°F	15 minutes	30 minutes <sup>1</sup>
Pre-vacuum	132°C / 270°F	4 minutes	

<sup>1</sup> Drying times vary according to load size and should be increased for larger loads.

### 5. Instructions for Use:

- A. Using a purpose designed drill or wire driver, drill the guide pin into the desired position.
- B. Place cannulated drill or reamer with an appropriately sized cannulation over guide pin.
- C. Drill or ream over guide pin ensuring to maintain a coaxial orientation with guide pin.
- D. Remove drill or reamer from pin.
- E. If using to pass suture; drill through bone, pass suture through eyelet, pull pin with attached suture through bone.
- F. Remove and dispose of pin.

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