## Parcus Nitinol Guide Wires with and without Needle Introducer

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## Parcus Nitinol Guide Wires with and without Needle Introducer (English)

### 1. Indications:

Parcus Nitinol Guide Wires with and without Needle Introducer are used to aid in the placement of cannulated implants and instruments.

### 2. Warnings:

- A. Caution: Federal Law restricts this device to sale by or on the order of a physician.
- B. This product is designed for Single Use Only.
- C. Inspect for damage prior to use. Do not use if damaged.
- D. Devices that are reused can fatigue and break.
- E. Do not use if the packaging has been opened or compromised.
- F. This product is not intended for use as an implant.
- G. This device is sold as Non-sterile or Sterile. Non-sterile devices must be sterilized prior to use. (See Sterilization)

### 3. Material:

Guide Wires are manufactured from Nitinol per ASTM 2063.

Needles and Stylets are made from stainless steel alloy, nickel-plated brass and plastic.

Materials used in the manufacture of this device that are intended to be placed inside the body are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

### 4. Sterilization:

- A. Nitinol Guide Wires with and without Needle Introducer may be supplied as Non-sterile and Sterile.

- B. Non-sterile devices:

The following table provides recommended minimum sterilization parameters that have been validated by Parcus Medical to provide a  $10^{-6}$  sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Pre-vacuum	132°C / 270°F	4 minutes	30 minutes <sup>1</sup>

<sup>1</sup> Drying times vary according to load size and should be increased for larger loads.

- C. Sterile Devices:

The device is supplied in sterile packaging. The contents are sterilized by EO gas. These products must never be re-sterilized.

### 5. Packaging and Labeling:

- A. Do not use this product if the packaging or labeling has been damaged, shows signs of exposure to moisture or extreme temperature or has been altered in any way.
- B. Please contact Parcus Medical Customer Service to report any package damage or alterations.

### 6. Storage:

Products must be stored in the original, unopened package in a dry place and must not be used beyond the expiration date indicated on the package.

### 7. Instructions for use:

For Percutaneous Procedures:

- A. Introduce the Needle with Stylet through the skin and underlying tissue into the desired position.
- B. Remove the Stylet, insert the Nitinol Guide Wire and remove the Needle while leaving the Nitinol Guide Wire in position.
- C. Slide the cannulated implant or instrument over the Nitinol Guide Wire into the desired position.
- D. Remove the Nitinol Guide Wire and discard the Nitinol Guide Wire, Needle and Stylet.

For Open Procedures:

- A. Place the Nitinol Guide Wire in the desired position.
- B. Slide the cannulated implant or instrument over the Nitinol Guide Wire into the desired position.
- C. Remove and discard Nitinol Guide Wire.

**Note:** Open procedure instructions may vary due to instrument use, implant selection and application. This example illustrates a typical surgical scenario and is provided to show format, not content.

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