Parcus Threaded Dilator

Important Product Information

Produktinformationsblatt

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Parcus Threaded Dilator (English)

1. Indications:

The Parcus Threaded Dilator is designed to create a starter hole to ease the insertion of Interference Screws in Bone-Tendon/Ligament-Bone ACL Reconstructions. The Threaded Dilator is cannulated and is compatible with 1.5mm or less Guide Wires.

2. Warnings:

- A. Caution: Federal Law restricts this device to sale by or on the order of a physician.
- B. Do not use this product if it is bent or there is visible damage to the shaft or drive section.
- C. This device is only to be used when sufficient bone quality and implant/graft size exists such that fixation is not compromised.
- D. This device is not intended to be used as an implant.

3. Material:

This device is manufactured from a Stainless Steel Alloy and may include an aluminum alloy handle. The materials used in the manufacture of this device are radio-opaque and can, therefore, be detected with conventional X-Ray or fluoroscopy.

4. Inspection:

- A. Inspect the device for damage at all stages of handling.
- B. If damage is detected, consult the manufacturer for guidance.

5. Cleaning:

- A. Immediate rinsing and cleaning after use with an enzymatic detergent will effectively remove and prevent drying of adherent blood, tissue, etc.
- B. Scrub device with a soft brush, paying special attention to areas where debris might accumulate. Always avoid any harsh materials that can scratch or mar the surface of the device.
- C. Rinse the device thoroughly with water following the cleaning process.

6. Sterilization:

This device is not provided sterile. The following table provides recommended minimum sterilization parameters that have been validated by Parcus Medical to provide a 10⁻⁶ sterility assurance level (SAL):

Cycle Type	Minimum Temperature	Minimum Exposure Time (Wrapped)	Minimum Dry Time
Pre-vacuum	132°C / 270°F	4 minutes	30 minutes ¹

¹ Drying times vary according to load size and should be increased for larger loads.

7. Instructions for use:

- A. Place the Threaded Dilator over an appropriately placed 1.5mm Guide Wire to guide the dilator into the insertion site.
- B. Rotate the Threaded Dilator clockwise until the dilator has reached the desired depth.
- C. Rotate the Threaded Dilator counterclockwise to remove the dilator.
- D. Place the desired Interference Screw over the Guide Wire and insert per applicable Instructions for Use.

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