

Knotless PEEK CF Push-In Suture Anchors

Surgical Innovation | Customer Driven

Knotless PEEK CF Push-In Suture Anchors provide a step-saving alternative to conventional "knotted" suture anchors and eliminate the "knot stacks" that can be associated with soft tissue irritation or impingement.

Made of Carbon Fiber-reinforced PEEK-OPTIMATM, Knotless PEEK CF Push-In Suture Anchors are non-absorbable, radiolucent, and MR safe with a modulus of elasticity closely matching cortical bone. Recommended for use in both large- and small-joint repairs.

Key Features

- · Diameter options: 2.8mm, 3.5mm, 4.5mm, & 5.5mm
- · Single-use and reusable drivers available
- · Suture passer included
- Made of Carbon Fiber-reinforced PEEK-OPTIMA from Invibio

Knotless PEEK CF Push-In Suture Anchors		
Part #	Description	
10312	28 Knotless PEEK CF Push-In Suture Anchor w/Driver, w/Suture Passer	
10994	35 Knotless PEEK CF Push-In Suture Anchor w/Driver, w/Suture Passer	
10313	35 Knotless PEEK CF Push-In Suture Anchor w/Suture Passer	
10314	45 Knotless PEEK CF Push-In Suture Anchor w/Suture Passer	
10315	55 Knotless PEEK CF Push-In Suture Anchor w/Suture Passer	

Required Instrumentation for Knotless PEEK CF Push-In Suture Anchors		
Part #	Description	
10344	35 Knotless Suture Anchor Driver (Use w/10313 Suture Anchor)	
10483	35 Knotless Suture Anchor Driver, Unsterile Reusable (Use w/10313 Suture Anchor)	
10345	45 or 55 Knotless Suture Anchor Driver (Use w/10314 or 10315 Suture Anchor)	
10487	45 or 55 Knotless Suture Anchor Driver, Unsterile Reusable (Use w/10314 or 10315 Suture Anchor)	
10365	3.2mm Awl w/Handle	
10538U	2.8mm x 300mm Drill Bit w/Positive Stop, Unsterile Reusable	
10615U	28 Knotless Drill Bit, 300mm, w/Positive Stop, Unsterile Reusable	
10537U	3.0mm x 300mm Drill Bit w/Positive Stop, Unsterile Reusable	
10823	8-Point Drill Guide w/Beaver Tail Handle, 100mm Length	
10924	8-Point Slotted Knotless Guide w/Beaver Tail Handle, Accommodates 300mm Drill Guide	
10446	8-Point Drill Guide w/Obturator, 300mm Length	
10330	V-Mouth Drill Guide w/Obturator, 300mm Length	
10406	Blunt Tip Obturator (Use w/10446 & 10330 Drill Guides)	
10771	Cannulated Conical Tip Obturator, 1.8mm Cannulation (Use w/10330 or 10446 Drill Guide)	
10405	Trocar Tip Obturator (Use w/10446 & 10330 Drill Guides)	
10397	Switching Stick	





PHONE 1.941.755.7965 • FAX 1.941.755.6543 • 6423 PARKLAND DRIVE SARASOTA, FL 34243 • www.parcusmedical.com