



V-LoX and V-LoX³ PEEK CF Screw-In Suture Anchors, the first family of suture anchors manufactured from carbon fiber-reinforced PEEK-OPTIMA, are recommended for use in both large and small-joint repairs.

Features & Benefits

Fully-threaded design

- Optimal cortical fixation
- Reduce the risk of "pull-back"

Made of carbon fiber-reinforced PEEK-OPTIMA

- Non-absorbable, radiolucent, and MR safe
- Modulus of elasticity closely matching cortical bone

V-LoX PEEK CF Screw-In Suture Anchors

Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
10253	5.5	12.7	V-LoX, screw-in, suture anchor, w/2, #2 sutures (wht/blu & wht/blk)	PEEK CF	1
10253T	5.5	12.7	V-LoX, screw-in, suture anchor, w/2, #2 sutures (blu & wht/grn)	PEEK CF	1
10254T	6.5	12.7	V-LoX, screw-in, suture anchor, w/2, #2 sutures (blu & wht/grn)	PEEK CF	1
10263	5.5	12.7	V-LoX, screw-in, suture anchor, w/2, #2 sutures (wht & wht/blu), w/needles (MO-6)	PEEK CF	1
10264	6.5	12.7	V-LoX, screw-in, suture anchor, w/2, #2 sutures (wht & wht/blu), w/needles (MO-6)	PEEK CF	1

V-LoX³ PEEK CF Screw-In Suture Anchors

Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
10324T	5.5	12.7	V-LoX ³ , screw-in, suture anchor, w/3, #2 sutures (blu, wht/grn & wht/blk)	PEEK CF	1
10325T	6.5	12.7	V-LoX ³ , screw-in, suture anchor, w/3, #2 sutures (blu, wht/grn & wht/blk)	PEEK CF	1

V-LoX PEEK CF Screw-In Suture Anchors w/Suture Tape

Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
11026	5.5	12.7	V-LoX, screw-in, suture anchor, w/2, 1.6mm suture tapes (wht/blu & wh/blk)	PEEK CF	1
11027	6.5	12.7	V-LoX, screw-in, suture anchor, w/2, 1.6mm suture tapes (wht/blu & wh/blk)	PEEK CF	1

Instrumentation for V-LoX & V-LoX³ PEEK CF Screw-In Suture Anchors

Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
11000U	4	Small Joint	Drill bit	Unsterile	Reusable
11241	4	Small Joint	Drill guide	Unsterile	Reusable
10338	4.1	Shoulder	Awl	Unsterile	Reusable
10277	5.5	Shoulder	5.5mm or 6.5mm punch/tap	Unsterile	Reusable
10222	5.5	Shoulder	5.5mm or 6.5mm punch/tap, hudson	Unsterile	Reusable



Shoulder