



MiT<sub>i</sub> and MiTi XL Screw-In Suture Anchors are recommended to address a variety of small-joint applications due to their small diameter and tremendous pull-out strength. The MiTi XL, featuring an extended shaft, is ideal for the arthroscopic management of shoulder instability.

## Features & Benefits

### Fully-threaded design

- Optimal cortical fixation
- Reduce the risk of "pull-back"

### Made of 6Al-4V ELI titanium alloy

- Self-tapping threads

MiTi Titanium Screw-In Suture Anchors					
Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
10425	2	6.4	MiT <sub>i</sub> , screw-in, suture anchor, w/2, 3-0 sutures (wht & wht/blu), w/needles (BB-1), w/drill bit	Titanium	1
10433	2.5	6.4	MiT <sub>i</sub> , screw-in, suture anchor, w/1, #2 suture (wht/blu), w/needles (MO-6), w/drill bit	Titanium	1
10429	2.5	6.4	MiT <sub>i</sub> , screw-in, suture anchor, w/2, 3-0 sutures (wht & wht/blu), w/needles (BB-1), w/drill bit	Titanium	1
10472	3.5	7.6	MiT <sub>i</sub> , screw-in, suture anchor, w/1, #2 suture (wht/blu), w/needles (MO-6), w/drill bit	Titanium	1

MiTi XL Arthroscopic Length Titanium Screw-In Suture Anchor					
Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
10434	2.5	6.4	MiT <sub>i</sub> XL, screw-in, suture anchor, w/1, #2 suture (wht/blu), w/arthroscopic length driver	Titanium	1

Instrumentation for MiTi XL Arthroscopic Length Titanium Screw-In Suture Anchors					
Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10446	-	Shoulder	8-Point drill guide, w/trocar tip obturator	Unsterile	Reusable
10330	-	Shoulder	V-Mouth drill guide w/trocar tip obturator	Unsterile	Reusable
10406	3.4	Shoulder	Blunt tip obturator <sup>1</sup>	Unsterile	Reusable
10771	3.3	Shoulder	Conical tip obturator, w/1.8mm cannulation	Unsterile	Reusable
10405	3.4	Shoulder	Trocar tip obturator <sup>1</sup>	Unsterile	Reusable



**Shoulder**