



Parcus Braid, available in multiple configurations: suture tape, infinity loop, and suture, is recommended for use in soft tissue approximation and allograft tissues in sports medicine procedures.

Features & Benefits

Made of Ultra High Molecular Weight Polyethylene (UHMWPE)

- High strength
- Non-absorbable

Variety of color combinations

- Suture management in complicated repairs

Various needle options

- Ideal for open surgical procedures

Parcus Braid Suture					
Part #	Size	Description	Total Length (in)	Working/Loop Length (in)	QTY/ Box
20046	#1	Parcus Braid, suture (wht/blu)	36	-	12
20050	#2	Parcus Braid, suture (wht/blu)	36	-	12
20049	#2	2, Parcus Braid, sutures (wht/blu & wht/blk)	36	-	12
20049S	#2	2, Parcus Braid, sutures (wht/blu & wht/blk)	36	-	1
20035S	#2	Parcus Braid, suture (blu), w/needles (MO-6 & CE)	36	-	1
20055	#2	Parcus Braid, suture (wht), w/needles (MO-6)	36	-	12
20055S	#2	Parcus Braid, suture (wht), w/needles (MO-6)	36	-	1
20054	#2	Parcus Braid, suture (wht/blk), w/needles (MO-6)	36	-	12
20054S	#2	Parcus Braid, suture (wht/blk), w/needles (MO-6)	36	-	1
20059	#5	Parcus Braid, suture (wht), w/needles (CCS)	36	-	12
20059S	#5	Parcus Braid, suture (wht), w/needles (CCS)	36	-	1

Parcus Braid Infinity Loop					
Part #	Size	Description	Total Length (in)	Working/Loop Length (in)	QTY/ Box
20053	#2	Parcus Braid, infinity loop (wht), w/needle	24	20	12
20053S	#2	Parcus Braid, infinity loop (wht), w/needle	24	20	1
20056	#2	Parcus Braid, infinity loop (wht/blu), w/needle	24	20	12
20056S	#2	Parcus Braid, infinity loop (wht/blu), w/needle	24	20	1
20060	#5	Parcus Braid, infinity loop (wht/blu), w/needle	24	20	12

Parcus Braid Suture Tape					
Part #	Size	Description	Total Length (in)	Working/Loop Length (in)	QTY/ Box
20032S	1.6mm	Parcus Braid, suture tape (wht/blu)	39	27	1
20034S	1.6mm	Parcus Braid, suture tape (wht/blk)	39	27	1
20042	2mm	Parcus Braid, suture tape (wht/blu)	39	7	6
20042S	2mm	Parcus Braid, suture tape (wht/blu)	39	7	1



Key: Needles	
Type	Description
CE	Reverse cutting curved needle
MO-6	Curved, 28mm needle
MO-6 & CE	MO-6 needle on one side, CE needle on the other side
BB-1	Curved 17mm needle with taper point
CCS	48mm curved piercing needle