

SLiK Fix is a screw-in tenodesis system recommended for use in both large and small-joint repairs. **SLiK Fix** is ideally suited for biceps repair. Available in both carbon fiber-reinforced PEEK-OPTIMA and natural PEEK-OPTIMA.

Features & Benefits

Activation button

- Allows for one-handed insertion which saves time in the operating room

Suture retention tip and tendon washer

- Effective manipulation of tendon

Designed to advance tendon prior to fixation


- Compression fixation along full length of implant

SLiK Fix PEEK CF Screws

Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
10995	7	15	SLiK Fix , tenodesis screw	PEEK CF	1
10996	8	15	SLiK Fix , tenodesis screw	PEEK CF	1
10997	9	15	SLiK Fix , tenodesis screw	PEEK CF	1

Instrumentation for SLiK Fix Screws

Part #	Diameter (mm)	Joint	Description	Sterile/ Unsterile	Single-use/Reusable
11016U	7	Shoulder	SLiK Fix , drill bit	Unsterile	Single-use
11017U	8	Shoulder	SLiK Fix , drill bit	Unsterile	Single-use
11018U	9	Shoulder	SLiK Fix , drill bit	Unsterile	Single-use
11019U	10	Shoulder	SLiK Fix , drill bit	Unsterile	Single-use



Shoulder