

PEEK CF Push-In Suture Anchors are recommended for knotted fixation in shoulder instability or foot and ankle.

Features & Benefits

Made of carbon fiber-reinforced PEEK-OPTIMA


- Non-absorbable, radiolucent, and MR safe
- Modulus of elasticity closely matching cortical bone

PEEK CF Push-In Suture Anchors

Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
10323T	3.5	8	Push-in, suture anchor, w/ 2, #2 sutures (blu & wht/grn)	PEEK CF	1
10407	3.5	8	Push-in, suture anchor, w/ 2, #2 sutures (blu & wht/grn), w/needles (MO-6 short)	PEEK CF	1

Instrumentation for PEEK CF Push-In Suture Anchors

Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10333U	2.8	Small Joint	Drill bit, w/positive stop	Unsterile	Reusable
10538U	2.8	Shoulder	Drill bit w/positive stop	Unsterile	Reusable
10482U	3	Small Joint	Drill bit w/positive stop	Unsterile	Reusable
10537U	3	Shoulder	Drill bit w/positive stop	Unsterile	Reusable
10823	-	Small Joint	8-Point drill guide, w/beaver tail handle	Unsterile	Reusable
10924	5.2	Shoulder	8-Point slotted knotless drill guide, w/beaver tail handle ¹	Unsterile	Reusable
10446	-	Shoulder	8-Point drill guide w/trocar tip obturator	Unsterile	Reusable
10330	-	Shoulder	V-Mouth drill guide w/trocar tip obturator	Unsterile	Reusable
10406	3.4	Shoulder	Blunt tip obturator ²	Unsterile	Reusable
10771	3.3	Shoulder	Conical tip obturator, w/1.8mm cannulation ²	Unsterile	Reusable
10405	3.4	Shoulder	Trocar tip obturator ²	Unsterile	Reusable
10397	2.8	Shoulder	Switching stick	Unsterile	Reusable



Shoulder

¹Accommodates 300mm drill guide ²Use w/10446 & 10330