

Draw Tight Suture-Based Anchors are recommended for use in both large and small-joint repairs. Made with UHMWPE and a natural PEEK-OPTIMA power tip.

## Features & Benefits

### Small insertion footprint & sub-cortical fixation

- Bone sparing

### Deployment sutures

- Tactile confirmation the anchor has been successfully deployed

### PEEK power tip

- Delivers improved insertion and a solid core for stronger fixation

#### Draw Tight DEX Suture-Based Anchors

Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
11034	1.8	-	Draw Tight DEX, push-in, suture anchor, w/1, 1.6mm suture tape (wht/blu), w/needles (MO-6)	UHMWPE	1

#### Draw Tight Suture-Based Anchors


Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
10565	1.8	-	Draw Tight, push-in, suture-based anchor, w/1, #2 suture (wht/blu)	UHMWPE	1
10820	3.2	-	Draw Tight, push-in, suture-based anchor, w/2, #2 sutures (blu & blu/blu)	UHMWPE	1
10588	3.2	-	Draw Tight, push-in, suture-based anchor, w/2, #2 sutures (wht/blu & wht/blk)	UHMWPE	1

#### Draw Tight Suture-Based Anchors w/Suture Tape

Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
11148	1.8	-	Draw Tight, push-in, suture-based anchor, w/1 1.6mm suture tape (wht/blu)	UHMWPE	1
11151	3.2	-	Draw Tight, push-in, suture-based anchor, w/2, 1.6mm suture tapes (wht/blu & wht/blk)	UHMWPE	1
11288	3.2	-	Draw Tight, push-in, suture-based anchor, w/2, 1.6mm suture tapes (blu, wht/blu)	UHMWPE	1
11289	3.2	-	Draw Tight, push-in, suture-based anchor, w/2, 1.6mm suture tapes (blk, wht/blk)	UHMWPE	1

#### Instrumentation for Draw Tight Suture-Based Anchors

Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10365	3.2	Shoulder	Awl	Unsterile	Reusable
10873U	-	Small Joint	1.8 Draw Tight, drill bit, w/positive stop	Unsterile	Reusable
10590U	-	Shoulder	1.8 Draw Tight, drill bit, w/positive stop	Unsterile	Reusable
10505U	3.2	Small Joint	Drill bit, w/positive stop	Unsterile	Reusable
10540U	3.2	Shoulder	Drill bit, w/positive stop	Unsterile	Reusable
10591	-	Shoulder	3.2 Draw Tight, drill guide <sup>1</sup>	Unsterile	Reusable
10823	-	Small Joint	8-point drill guide, w/beaver tail handle	Unsterile	Reusable
10446	-	Shoulder	8-point drill guide, w/trocar tip obturator	Unsterile	Reusable
10406	3.4	Shoulder	Blunt tip obturator <sup>2</sup>	Unsterile	Reusable
10771	3.3	Shoulder	Conical tip obturator, 1.8mm cannulation <sup>2</sup>	Unsterile	Reusable
10405	3.4	Shoulder	Trocar tip obturator <sup>2</sup>	Unsterile	Reusable



Shoulder

<sup>1</sup>Use w/10329 or 10466    <sup>2</sup>Use w/10446 & 10330