

Knotless PEEK Push-In Suture Anchors are recommended for use in both large and small-joint repairs. Available in natural PEEK-OPTIMA.


Features & Benefits

Knotless

- Provides a step-saving alternative to conventional "knotted" suture anchors
- Eliminates "knot stacks" associated with soft tissue irritation

Made of PEEK-OPTIMA Natural

- Non-absorbable, radiolucent, MR safe



Knotless PEEK Push-In Suture Anchors

Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
11321	2.8	10	28 Knotless, push-in, suture anchor, w/suture passer	PEEK	1

Knotless PEEK Push-In Suture Anchor, Hip Length

Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
11321H	2.8	10.2	28 Knotless, push-in, suture anchor, w/hip length driver, w/suture passer	PEEK	1

Instrumentation for Knotless PEEK Push-In Suture Anchors

Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
10874U	-	Small Joint	28 Knotless drill bit, w/positive stop	Unsterile	Reusable
10615U	-	Shoulder	28 Knotless drill bit, w/positive stop	Unsterile	Reusable
10823	-	Small Joint	8-point drill guide, w/beaver tail handle	Unsterile	Reusable
10924	5.2	Shoulder	8-point slotted knotless drill guide, w/beaver tail handle	Unsterile	Reusable
10446	-	Shoulder	8-point drill guide, w/trocar tip obturator	Unsterile	Reusable
10330	-	Shoulder	V-Mouth drill guide w/trocar tip obturator	Unsterile	Reusable
10406	3.4	Shoulder	Blunt tip obturator	Unsterile	Reusable
10771	3.3	Shoulder	Conical tip obturator, w/1.8mm cannulation	Unsterile	Reusable
10405	3.5	Shoulder	Trocar tip obturator	Unsterile	Reusable



Shoulder