



Knotless AP Push-In Suture Anchors are recommended for use in both large and small-joint repairs. The proven design of the Knotless Product family, is updated with Parcus Advanced Polymer (AP) material that is comprised of 70% PLGA and 30% B-TCP.

Features & Benefits

Advanced Polymer

- Resorbable, radiolucent, and MR safe

Knotless

- Provides a step-saving alternative to conventional "knotted" suture anchors
- Eliminates "knot stacks" associated with soft tissue irritation

Knotless AP Push-In Suture Anchors

Part #	Diameter (mm)	Length (mm)	Description	Material	QTY/Box
10802	2.8	10.1	28 Knotless, push-in, suture anchor, w/suture passer	AP	1
10803	4.5	16.3	45 Knotless, push-in, suture anchor, w/suture passer	AP	1
10804	5.5	18.1	55 Knotless, push-in, suture anchor, w/suture passer	AP	1

Instrumentation for Knotless AP Push-In Suture Anchors

Part #	Diameter (mm)	Joint	Description	Sterile/Unsterile	Single-use/Reusable
11000U	4	Small Joint	Drill bit	Unsterile	Reusable
11241	4	Small Joint	Drill guide	Unsterile	Reusable
10365	3.2	Shoulder	Awl	Unsterile	Reusable
10874U	-	Small Joint	28 Knotless drill bit, w/positive stop	Unsterile	Reusable
10615U	-	Shoulder	28 Knotless drill bit, w/positive stop	Unsterile	Reusable
10823	-	Small Joint	8-point drill guide, w/beaver tail handle	Unsterile	Reusable
10924	5.2	Shoulder	8-point slotted knotless drill guide, w/beaver tail handle	Unsterile	Reusable
10446	-	Shoulder	8-point drill guide, w/trocar tip obturator	Unsterile	Reusable
10330	-	Shoulder	V-Mouth drill guide w/trocar tip obturator	Unsterile	Reusable
10406	3.4	Shoulder	Blunt tip obturator	Unsterile	Reusable
10771	3.3	Shoulder	Conical tip obturator, w/1.8mm cannulation	Unsterile	Reusable
10405	3.5	Shoulder	Trocar tip obturator	Unsterile	Reusable
11172	-	Shoulder	Awl, w/ taper from 3.6mm to 4.1mm	Unsterile	Reusable



Shoulder